

Introducing a hot or cold therapy eye mask to soothe and bring relief from:

- Dry Eye Symptoms
- Headache
- Sinus Pressure
- Tension
- Swelling

Oasis

REST & RELIEF™

EYE MASK

Ref# OM8000

Elastic strap for a comfortable fit

Adjustable closure for compression

Durable temperature-retaining silica beads for hot or cold therapy

Soft and reusable with contouring-fit for full contact and comfort



OASIS®  True Dry Eye Solutions.

Call to Order **(844) 820-8940**

Register: www.oasismedical.com

Email: customerservice@oasismedical.com

From the company that
brought your practice

Oasis
TEARS®

Includes washable cover.

- Soft, breathable fabric.



INSTRUCTIONS FOR USE:

INTENDED USE:

Hot/Cold therapy eye mask to soothe and bring relief from:

- Dry Eyes • Headache • Sinus Pressure • Tension • Swelling

HOT THERAPY

Directions:

1. Eye mask must be at room temperature before heating
2. Place eye mask flat on clean microwave safe dish.
3. Set microwave power to high and heat for 20 seconds. Before applying, always touch eye mask with fingers first to test for desired temperature. If eye mask is too hot, allow it to cool for 2 minutes and test again. For additional heat, add in 5 second increments to obtain your desired temperature. Do not exceed a maximum of 30 seconds of heating.
4. Place eye mask over closed eyes. Optionally secure adjustable strap.
5. Wear for 3 to 5 minutes or as directed by a physician. Do not apply for more than 20 minutes at a time.

COLD THERAPY

Directions:

1. Place eye mask on clean dry surface and place in the freezer for 30 minutes.
2. Remove mask from freezer and place over closed eyes. Wear for 5 to 7 minutes or as directed by a physician.

Care: Surface clean mask with a damp cloth. Do not submerge mask in water. Machine or hand wash cover only.

Storage: Store in cool dry place. Keep away from flames.

QUESTIONS

(844) 820-8940 or (909) 305-5400
Visit: www.oasismedical.com

WARNINGS AND PRECAUTIONS

Please read all information before use.

- For external use only.
- Keep out of reach of children.
- Do not apply for more than 20 minutes at a time.
- Dispose of mask if punctured or broken.
- Test eye mask temperature prior to use.
- Overheating or prolonged application may cause injury.
- Do not use on infants.
- Do not apply to sensitive skin, skin with no sensation, or open wounds.
- Do not boil.
- Overheating in a microwave may damage eye mask.
- Do not ingest contents.
- Warm according to instructions. Do not apply if too hot.
- Do not fall asleep using eye mask.
- People with diabetes, circulation disorders, skin infections, etc. should consult with their doctor before use.
- Consult your doctor if your condition worsens or irritation develops.

Other products by OASIS®:

Preservative-Free Lubricant Eye Drops:

- Oasis TEARS® - Single use vials
- Oasis TEARS® PLUS - Single use vials
- Oasis TEARS® PF - Bottle

Hydrating Eyelid & Lash Cleaners for Daily Hygiene:

- Oasis LID & LASH® + Tea Tree - Oil Wipes

