ENGLISH

GENERAL DESCRIPTION	2
UNPACKING INSTRUCTIONS.....	2
BATTERY SPECIFICATION AND INSTALLATION.....	3
INSTRUCTIONS FOR USE.....	4
SPECIFICATIONS	5
WARRANTY	6
SERVICE AND REPAIR	6
SUGGESTED OPERATION	7

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I. GENERAL DESCRIPTION

The Colvard Pupillometer is a precision instrument designed to measure pupillary dilation in low levels of illumination.

The light weight, hand held instrument is designed to fit comfortably in the user's hand, facilitating fast and accurate measurements. The Colvard Pupillometer is used by the examiner very much like a direct ophthalmoscope.

The Colvard Pupillometer has features designed to make the device safer, to protect the electronics of the instrument, and to extend the life of the lithium battery energy source.

The Colvard Pupillometer is:

1. **Easy to use.**

2. **Portable.**

An easy to store carry case is included with each Pupillometer.

3. **Battery operated.**

A lithium energy source operates the instrument.

Caution: This product contains Dry Natural Rubber which may cause allergic reactions.

II. UNPACKING INSTRUCTIONS

- A. Remove all protective packaging materials.
- B. Check for missing items. The following items should be present:
 - Protective foam-lined case
 - Colvard Pupillometer
 - Fresh 3 volt lithium battery
 - Instruction Manual
- C. Visually check the instrument and accessories for damage.
- D. Notify OASIS Medical Customer Service immediately if any components are missing or damaged.

Figure 1

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III. BATTERY SPECIFICATION AND INSTALLATION

A. Battery Specification:

The Pupillometer uses one 3.0 volt lithium battery (type CR123A). Replacement batteries can be obtained from a photographic dealer or ordered from OASIS Medical.

NOTE: The Pupillometer will not function properly with mercury batteries.

B. Battery Installation:

1. Locate the battery compartment (see Figure 2). Open the battery compartment by twisting off the cap.
2. Insert the lithium battery into the compartment according to the diagram on the pupillometer handle.
3. Replace the cap.

NOTE: If the pupillometer does not work, check that the battery was inserted correctly.

CAUTION: Do not disassemble, crush, recharge, or expose batteries to fire. Dispose in accordance with appropriate regulations. Opened batteries should be treated as a hazardous waste.


Figure 2. Bottom view - battery compartment

CAUTION

CHECK TO BE SURE THAT THE BATTERY IS INSTALLED CORRECTLY. INCORRECT INSTALLATION COULD CAUSE SEVERE DAMAGE TO THE ELECTRONICS AND INVALIDATE THE WARRANTY.

IV. INSTRUCTIONS FOR USE

A. Exam room preparation:

To approximate the level of illumination encountered during night driving, turn off all the lights in the examination room and “crack” the hall door slightly (1-2 inches). There should be enough light to allow the technician to examine the patient safely.

B. Operation:

The Colvard Pupillometer is activated by pressing the ON/OFF power switch located on the instrument handle (see Figure 3).

As long as the ON/OFF button is held down, the bright enhanced image continues. When switch is released, the enhanced image will persist briefly, but the quality of the image will diminish rapidly.


Figure 3. Side view

C. Pupil measurement:

During examination, instruct the patient to focus on an object at the end of the examination room with the fellow eye. The examiner should position himself/herself so as not to obstruct the patient line of sight with the fellow eye.

Look through the eyepiece to visualize the iris and pupil. Move the instrument slightly forward or back to bring the image into sharper focus. The enhanced image is normally somewhat grainy with a green phosphorescent glow. The image should be clear enough to visualize the pupil with ease.

A reticule in the device superimposes a millimeter ruler over the image of the iris and pupil and allows for easy measurement of the pupillary diameter.

Record your measurement.

Repeat the procedure with the second eye.

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
Figura 3. Vista frontale

D. BEFORE turning on the exam room light:

The electronics of the device contains a protective mechanism which helps to shield the sensitive light amplification instrumentation from sudden exposure to high levels of light energy. However, care should be taken to release the light amplification activation button before the device is exposed to ordinary room illumination.

E. Recommended cleaning methods and materials:

- Clean only with a soft cloth. Cloth may be moistened with isopropyl alcohol. Dry thoroughly with a soft cloth.
- Clean lenses with commercially available lens cleaner only.

F. Storage:

Store the Colvard Pupillometer in its protective foam-lined case when not in use. Keep the instrument dry at all times. Moisture or excessive humidity could damage the electronics. Avoid temperatures below freezing and above 65°C (150°F).

V. SPECIFICATIONS

A. Classification: Type B

Type B medical electronic equipment provides protection against electric shock, particularly in regards to allowable Leakage Current.

The IR illumination diode is classified as a Class I LED Product.

B. Electrical/Mechanical:

Power Supply – 3.0-volt lithium battery (type CR123A)

Light Amplification – Image Intensifier Tube

Illumination IR Diode – Peak Wavelength: 810 nm

Focal Length of Front Optics: 32 mm

Reticle Geometry: Horizontal Axis – 13 mm 1 mm graduations; Vertical Axis – 13 mm 1 mm graduations from 0mm to 4mm and from 9mm to 13mm. (Implemented June 2002)

Accuracy of Reticle: +/- 0.1 mm

WARNINGS!! DO NOT OPEN INSTRUMENT CASE!!

This instrument contains a High Voltage power supply inside the protective case of the pupillometer which takes at least 10 minutes to fully discharge.

The image intensifier tube can be easily damaged by touching its surfaces.

Do not operate pupillometer in high ambient light levels as this will severely shorten the life of the Image Intensifier Tube.

VI. WARRANTY

OASIS Medical warrants its new equipment to be free from defects in workmanship or materials. Any product which is proven to be defective in workmanship or materials will be repaired or replaced at our discretion, free of charge, up to one year from date of purchase. This warranty covers all repairs and service of parts that have proved defective by manufacture and not by use or mishandling. This type of service will be handled in our home office. Shipping charges for returns for repair of non-warranted items will be the responsibility of the customer. Alteration, repair or modification of any product which is performed by persons not authorized by OASIS Medical, Inc. will result in immediate loss of warranty.

VII. SERVICE AND REPAIR

Technical information and operational instruction may be obtained by calling OASIS Medical Customer Service at (800) 528-9786 or fax (909) 305-9987.

When it is deemed necessary to return your Pupillometer for factory service, you must obtain a Return Goods Authorization (RGA) number from OASIS Customer Service at the above number.

Items sent to OASIS for repair must have the issued RGA number printed on the outside of the package. The package should be addressed as follows:

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SUGGESTED OPERATION OF THE COLVARD PUPILLOMETER

As pupil size continues to play a significant role in the outcomes of refractive surgery, the pre-operative use of the Pupillometer continues to be the standard of care for obtaining accurate measurement. In an ongoing effort to standardize methods of pre-operative pupil measurement OASIS Medical has developed the following checklist that should be followed for the greatest level of accuracy and consistency.

- Obtain screening room light level of approximately 4 cd/m². It should be noted that the light level might need to be made darker depending on the lifestyle of the particular patient.
- Measure the light level in the room with a light meter
- Wait at least 30 seconds until taking the measurement for full dilation of the pupil
- Once the Pupillometer is placed over the eye, instruct the patient to focus on an object at the end of the room with the fellow eye.
- Focus the reticule and take the measurement
- Document measurement in patient records
- Repeat process on the second eye

Using the Colvard Pupillometer and following this checklist will enable the measurement to obtain a highest level of accuracy and repeatability. It should also be noted that these are suggestions that may need to be modified based on the individual circumstances of an individual case.

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Glendora, CA 91741 U.S.A.
Donawa Consulting Srl