

OASIS   
SURGICAL PRODUCTS

# Premier Edge® Microsurgical Knives - Disposable

**Ergonomic Handle** Designed for comfort, safety, efficiency and ease of use to provide the most natural feel in the surgeon's hand and enhance overall control during surgery

**Blade** Consistently sharp with non-glare finish creates reproducible results

**Packaging** Sterile blister packs, 6 per box

**STAB**



**SLIT  
THIN LINE**



**CLEAR CORNEA  
THIN LINE**



**PARALLEL SIDE**



**TRAPEZOID**



**IMPLANT**



**SIDE PORT**



**ROUND TUNNEL**



**CRESCENT TUNNEL**



**MVR**






REF #	DESCRIPTION	WIDTH	ANGLE	BEVEL
<b>STAB KNIVES</b>				
PE3015	Stab	15°		Double
PE3015-3	Stab	15°	3mm depth	
PE3015-5	Stab	15°	5mm depth	
PE3022	Stab	22.5°		Double
PE3030	Stab	30°		Double
PE3045	Stab	45°		Double
<b>THIN LINE SLIT KNIVES</b>				
PE4820-TL	Thin Line Slit	2.0mm	Angled	Bevel Up
PE4822-TL	Thin Line Slit	2.2mm	Angled	Bevel Up
PE4824-TL	Thin Line Slit	2.4mm	Angled	Bevel Up
PE4826-TL	Thin Line Slit	2.65mm	Angled	Bevel Up
PE4827-TL	Thin Line Slit	2.75mm	Angled	Bevel Up
PE4828-TL	Thin Line Slit	2.85mm	Angled	Bevel Up
PE4830-TL	Thin Line Slit	3.0mm	Angled	Bevel Up
PE4832-TL	Thin Line Slit	3.2mm	Angled	Bevel Up
<b>THIN LINE CLEAR CORNEA KNIVES</b>				
PE4122-TL	Thin Line Clear Cornea	2.2mm	Angled	Double
PE4125-TL	Thin Line Clear Cornea	2.5mm	Angled	Double
PE4126-TL	Thin Line Clear Cornea	2.65mm	Angled	Double
PE4127-TL	Thin Line Clear Cornea	2.75mm	Angled	Double
PE4128-TL	Thin Line Clear Cornea	2.85mm	Angled	Double
PE4130-TL	Thin Line Clear Cornea	3.0mm	Angled	Double
PE4132-TL	Thin Line Clear Cornea	3.2mm	Angled	Double
<b>CLEAR CORNEA PARALLEL SIDE</b>				
PE3226	Parallel Side	2.65mm	Angled	Double
PE3227	Parallel Side	2.75mm	Angled	Double
PE3228	Parallel Side	2.85mm	Angled	Double
PE3230	Parallel Side	3.0mm	Angled	Double
<b>CLEAR CORNEA TRAPEZOID</b>				
PE3302	Trapezoid	2.5mm – 3.0mm	Angled	Bevel Up
<b>IMPLANT KNIVES</b>				
PE3432	Implant	3.2mm	Angled	Bevel Up
PE3435	Implant	3.5mm	Angled	Bevel Up
PE3440	Implant	4.0mm	Angled	Bevel Up
PE3452	Implant	5.2mm	Angled	Bevel Up
<b>SIDE PORT KNIVES</b>				
PE3010	Side Port	1.0mm	Angled	Double
PE3012	Side Port	1.2mm	Angled	Double
PE3014	Side Port	1.4mm	Angled	Double
<b>ROUND TUNNEL KNIVES</b>				
PE3622	Round Tunnel	2.2mm	Angled	Double
PE3630	Round Tunnel	3.0mm	Angled	Double
<b>CRESCENT TUNNEL KNIVES</b>				
PE3720	Crescent Tunnel	2.0mm	Straight	Bevel Up
PE3720A	Crescent Tunnel	2.0mm	Angled	Bevel Up
PE3720D	Crescent Tunnel	2.0mm	Angled	Bevel Down
<b>MVR KNIVES</b>				
PE3519	MVR	19g	Straight	Double
PE3519A	MVR	19g	Angled	Double
PE3520	MVR	20g	Straight	Double
PE3520A	MVR	20g	Angled	Double
PE3523	MVR	23g	Straight	Double
PE3523A	MVR	23g	Angled	Double
PE3525	MVR	25g	Straight	Double



# Needles, Cannulas, Cystotomes – Disposable




■ Available sterile 10 per box\*

Gauge	Anesthesia	Visco/Irrigation Cannulas	Capsulotomy Instruments
19g (1.10)	 4005 (1"/25mm) 4005TP (1"/25mm)	 4036 (1"/25mm)	
20g (0.90)		4001V-27 (0.875"/22mm) Stepped down to 27g (0.40) 10mm from tip	
21g (0.80)			
23g (0.60)	4004C (1.5"/38mm) Atkinson 4004G (0.875"/22mm) Atkinson 4004F (1.25"/32mm) Atkinson	4001B (0.875"/22mm) 10mm Bend 4001CF (0.875"/22mm) 4mm Bend 4001CTW (0.875"/22mm) 8mm Thin Wall	4017J (0.625"/16mm) Formed
24g (0.55)			
25g (0.50)	4004 (1.5"/38mm) 4004B (1.5"/38mm) Atkinson 4004D (0.875"/22mm) Atkinson 4004E (0.625"/16mm) Atkinson 4004M (1.25"/32mm) Atkinson	4001A (0.875"/22mm) 6mm Bend 4001AF (0.875"/22mm) 4mm Bend 4001ATW (0.875"/22mm) 6mm Thin Wall 4003 (1.125"/29mm) Curved	4017 (0.875"/22mm) Formed 4017A (0.875"/22mm) Straight 4017G (0.625"/16mm) Formed 4017IE (0.625"/16mm) Formed 4017L (0.625"/16mm) Reverse Formed 4017M (0.625"/16mm) Angled 4017S (0.625"/16mm) Formed 4017Z (0.625"/16mm) F-SL Curve
26g (0.45)			
27g (0.40)		4002 (0.875"/22mm) 6mm Bend 4002A (0.875"/22mm) 11mm Bend 4002B (0.875"/22mm) 8mm Bend 4002F (0.875"/22mm) 4mm Bend 4002L (0.875"/22mm) 9mm Bend	4017B (0.875"/22mm) Formed 4017IS (0.625"/16mm) Straight 4017IO (0.625"/16mm) Formed 4017N (0.625"/16mm) Angled 4017R (0.625"/16mm) Formed 4017Y (0.625"/16mm) F-SL Curve
30g (0.30)		4001 (0.875"/22mm) 6mm Bend 4001F (0.875"/22mm) 4mm Bend	4017I (0.625"/16mm) Formed 4017IR (0.625"/16mm) Reverse Formed

View Surgical Technique Videos at [www.oasismedical.com](http://www.oasismedical.com)

## Cannulas, Cortex Extractors, Capsule Polishers – Disposable

■ Available sterile 10 per box\*

Gauge	Hydrodissection/LASIK Cannulas	Cortex Extractors and Vectis	Capsule Polishers
			
<b>19g (1.10)</b>			
<b>20g (0.90)</b>			
<b>21g (0.80)</b>		4016 (0.875"/22mm) U-shaped	4027 (0.875"/22mm) 9mm Bend 4027B (0.875"/22mm) 6mm Bend
<b>23g (0.60)</b>	4123 (0.875"/22mm)	4007 Simcoe Double I and A 4008 (0.5"/13mm) Simcoe 4009 (0.875"/22mm) Simcoe 4016F (0.875"/22mm) U-shaped 4026 (0.5"/13mm) Charleux	4027A (0.875"/22mm) 9mm Bend 4028B (0.875"/22mm) Silicone Sleeve
<b>24g (0.55)</b>		4029A Vectis, three 0.25mm ports	
<b>25g (0.50)</b>	4016AH (0.875"/22mm) U-shaped 4036A (0.875"/22mm) 8mm Bend 4036C (0.875"/22mm) 4mm Bend 4036F (0.875"/22mm) Helsinki 4036M (0.875"/22mm) Sauter 4068 (0.875"/22mm) Kellan 4036L (0.875"/22mm) Lin LASIK 4036LB (0.875"/22mm) Buratto LASIK 4036LR (0.875"/22mm) Lindstrom LASIK 4225 (0.875"/22mm) LASIK, 6 ports	4016G (0.875"/22mm) J-shaped	4027E (1"/25mm) Olive Tip
<b>26g (0.45)</b>	4036LS (1"/25mm) Slade LASIK		
<b>27g (0.40)</b>	4016BH (0.875"/22mm) U-shaped 4036B (0.875"/22mm) 8mm Bend 4036D (0.875"/22mm) 4mm Bend 4036G (0.875"/22mm) Helsinki 4036N (0.875"/22mm) Sauter 4036J (0.75"/19mm) Chang 4036P (0.875"/22mm) Cleaving, Round 4036Q (0.875"/22mm) Cleaving, Flat	4016B (0.875"/22mm) U-shaped	4028 (0.875"/22mm) Terry
<b>30g (0.30)</b>			


\* Also available in the U.S. through the Bulk Program


## Visco Shield® HPMC Viscoelastic

- Available in 20,000 and 8,000 cps
- Stores at room temperature
- Packaged sterile 6 per box with 23g x 7/8 thin wall cannula

## Visco Shield® TOPICAL

- Designed for topical use to hydrate and protect the outer surface of the eye
- No preservatives
- Packaged sterile 6 per box with 23g x 7/8 thin wall cannula

REF #	DESCRIPTION	
<b>Viscoelastic – 2ml syringe</b>		
5122	Visco Shield Viscoelastic 20,000 cps	
51082	Visco Shield Viscoelastic 8,000 cps	

REF #	DESCRIPTION	
<b>Topical – 1.5ml syringe</b>		
5200	Visco Shield Topical 8,000 cps	


## Iris Retractors – Disposable

- To temporarily stretch and retain the iris during surgery
- Flexible blue polypropylene hooks
- Flat silicone stopper adjusts easily to secure hook
- 6 pouches per box, 5 hooks each

REF #	DESCRIPTION	
4465	Iris Retractors Sterile, 30 per box	

## Disposable Bimanual Handpiece Sets

- Decreases risk of contamination
- For use with any phaco machine
- **Purple** Aspiration Handle
- **Blue** Irrigation Handle
- Thin wall cannula design ensures maximum flow
- Only 2 sideport incisions necessary
- Sterile, 10 sets per box



REF #	GAUGE	ASPIRATION TIP	IRRIGATION TIP
1719	19	Closed end, rough, with round 0.3mm port	Closed end, rough, with dual oval 0.5mm ports
1720	20		
1721	21		
1719-1	19	Closed end, rough, with round 0.3mm port	Open end, smooth
1720-1	20		
1721-1	21		
1719-S	19	Closed end, rough, with round 0.3mm port	Closed end, rough, with dual oval 0.5mm ports
1720-S	20		
1721-S	21		

## Handles and Bimanual Cannulas

REF #	HANDLE WITHOUT TIP	QUANTITY
4000-AS	Aspiration Handle, Sterile	10 per box
4000-ASV	Aspiration Handle, Vented, Sterile	10 per box
4000-IR	Irrigation Handle, Sterile	10 per box
REF #	BIMANUAL CANNULAS	QUANTITY
4121	Bimanual Irrigation Cannula 21g (0.80) x 0.875" (22mm), Sterile 0.51mm open end irrigation port	10 per box
4122	Bimanual Aspiration Cannula 23g (0.60) x 0.875" (22mm), Sterile 0.4mm top port. Micro etched 2mm from end for capsule polishing	10 per box

## Accessories – Disposable

	REF #	DESCRIPTION	QUANTITY
<b>Corneal Marking Pad and Pen</b>			
	0441	Corneal Marking Pad, Pre-Inked Gentian Violet Sterile.	10 per box
	0442	Corneal Marking Pen, Fine Tip, Gentian Violet Sterile.	10 per box
<b>Disc Fluid Filter</b>			
	0601	For use between syringe and needle. Filters out particulate and biological contaminate matter larger than 0.22 µ. Bidirectional filter: female luer lock and male luer slip with screw-in lock to secure cannula in place. Sterile.	50 per box
<b>Eye Cautery</b>			
	3501-A	Fine tip, battery operated. Low temperature. Sterile. Range: 1100°F ± 50°F	10 per box
<b>SOFT CELL® PVA Foam</b>			
	0525	Surgical Spears Sterile. High density, enhanced fluid capacity, fiber-free	150 per box
	0522-5	LASIK Eye Drain (Chayet) Sterile, 10mm internal diameter	20 per box
	0523-5	LASIK Eye Drain (Chayet/Katzen) Sterile, 11.5mm internal diameter	20 per box
	0529	Wick + Wipe 1 sterile wipe per pouch	20 per box

European Authorized Representative  
Donawa Consulting Srl  
Piazza Albania 10  
00153 Rome, Italy

OASIS® 

Glendora, CA  
909.305.5400 ■ 909.305.9987 FAX  
international@oasismedical.com  
www.oasismedical.com



ISO 13485:2003  
CM004 4/09