

SOFT SHIELD® Collagen Shield by OASIS® Medical



## Corneal Comfort & Protection

NOW CE APPROVED



- Ocular surface protection following surgery, injection, traumatic and non-traumatic corneal conditions
- Comfort of a thin, clear, and pliable collagen barrier
- Hydrates to 70% water content
- Protects the surface of the eye
- Lubricates the ocular surface with collagen while it gradually dissolves

Available in Quick Shield, 12 hour, 24 hour, and 72 hour dissolution times

Diameter: 14.5mm    Compound base curve: 9mm

OASIS®  True Surgical Solutions.

Call to Order

US Toll-Free: (844) 820-8940

International: (909) 305-5400

Email: [customerservice@oasismedical.com](mailto:customerservice@oasismedical.com)

Visit: [www.oasismedical.com](http://www.oasismedical.com)