

# VISION CARE SOLUTIONS



Lubricate



Compress



Eyelid  
Hygiene



Nutrition



Punctal  
Occlusion



MY OASIS

## QUALITY ASSURANCE WORLDWIDE

**The OASIS® Medical Quality System is structured to meet the requirements of the Food and Drug Administration (FDA) Quality System Regulations. OASIS® Medical is registered with the FDA as a:**

- Medical Device Establishment - OASIS® Quality System complies with the Code of Federal Regulations (21 CFR §820)
- Drug Establishment - OASIS® Quality System complies with the Code of Federal Regulations (21 CFR §211)

**The OASIS® Medical Quality Management System is certified by DQS Medizinprodukte GmbH to EN ISO 13485:2016, ISO 15223-1:2021, Medical Device Single Audit Program (MDSAP) and meets the requirements of Annex II - Excluding Section 4, of Council Directive 93/42/EEC concerning medical devices for applying the CE marking to products:**

CE Marked Product Lines



Lacrimal Plug (SOFT PLUG® Silicone Punctum Plugs (Pen-Style Inserter is pending CE mark), SOFT PLUG® Extended Duration Plugs, SOFT PLUG® Extended Duration 180 Canalicular Plugs, FORM FIT® Hydrogel Canalicular Plug)

Disposable Microsurgical Instruments  
Microsurgical Knives (PremierEdge® Microsurgical Knives)  
HPMC Viscoelastics (VISCO SHIELD® HPMC Viscoelastic)  
HPMC Viscoelastics - Visco Shield Topical (VISCO SHIELD® HPMC Viscoelastic Topical)  
PVA Foam Products (SOFT CELL® PVA Foam Sponge Products)

Surgical Marking Instruments (Surgical Marking Pen)  
Iris Expander (OASIS® Iris Expander)  
Lacrimal Instruments  
Collagen Corneal Shield (SOFT SHIELD® Collagen Shields)  
Iris Retractor




## 3 WAYS TO ORDER:

1. **Call:** (844) 820-8940
2. **Email:** customerservice@oasismedical.com
3. **Register online:** www.oasismedical.com

### Ask your OASIS Representative about The OASIS ECP Channel

This channel is a dedicated educational resource for understanding dry eye, punctal plugs, MGD, diagnostics, and contains ECP pearls. Videos are presented by guest speakers who share their first-hand experience on what is working and what isn't working for the benefit of the ECP Community serving Dry Eye Patients.

#### Featured Videos




**Reimbursements**

Guest Presenter: Jason Flores, OD

Title: Solving for the Dry Eye Equat... 57:24

20 DAYS AGO

Jason Flores OD, Solving for the...




**Reimbursements**

Guest Presenter: John Lahr, OD

Title: Coding, Billing, and Reimbursement Updates (New Insights for... 1:04:01

2 MONTHS AGO

John Lahr OD, Coding, Billing, a...



**Punctal Plugs**


Guest Presenter: Paul Karpecki, OD, FAAO

Title: Punctal Occlusion - Expand Your Knowledge, Your Protocols, and You... 1:04:45

2 MONTHS AGO

Paul Karpecki OD FAAO, Punct...

### Punctal Plug Options



Reimbursements Info

Punctal Plug "How To"

Patient Education

#### Reimbursement Videos

- The ECP Workshop Video: John Lahr OD, Coding, Billing, and Reimbursement Updates for 2020. Presented by John Lahr OD FAAO (New Insights for 2021)

#### "How To" Videos

- The ECP Workshop Video: Paul Karpecki OD, FAAO, Punctal Occlusion - Expand Your Knowledge, Your Protocols, and Your Practice

#### Patient Education Videos

- What is Dry Eye?
- VLOG: Patient's Journey with Punctal Plugs

### View the OASIS Punctal Plug Resource Page

Scan the Code



Easily identify information that may be helpful with supporting the implementation of punctal occlusion in your practice today.



# RESOURCE CENTER

If your practice is seeking alternative ways to bring vision care solutions that can make a difference to your dry eye sufferers, then OASIS would like to introduce to your practice the MY OASIS Platform.


My OASIS platform was developed for practices to support the expansion of their services through continuity of care. This automated platform provides a customized storefront, the ability to bridge and engage with your patient community, all while bringing the convenience of at-home delivery on recurring or non-recurring basis. Let the platform tools be your website marketing. Let your practice's key person reach new members and build your vision care / dry eye community, prevent drop off, and bring in supplemental income. MY OASIS cyber platform diversifies your outreach, which broadens your growth potential. It supports communication, connection to nurture relationships, and the ability gain new actionable insights.

### FOR THE DOCTOR, MY OASIS IS:

- A convenient approach to customizing and delivering a daily eye care regimen with OASIS solutions.
- Access to eye care products that make a difference through a custom OASIS® store-front. Invite your members to easily place orders remotely with your unique code (recurring order option available).
- The ability to invite members to learn through the educational information offered covering lubricating eye drops, eye lid hygiene and maintenance, how to use them, and more.
- Access to tools that display your practice's growth, patient retention, and at-home care compliance.






### FOR THE CONSUMER/DRY EYE SUFFERER, MY OASIS IS:

- Well being because of access to your at-home care.
- Consumers with dry eye are seeking a vision care provider that is committed to their eye care. MY OASIS puts your practice in the palm of your patient's hand.
- Having a resource center where they can learn about dry eye and educate themselves on their treatment options.
- Convenience of ordering their doctor recommended solutions from home.
- The ease of setting their at-home care on a scheduled shipment.
- The ability to stay up to date on new products and services.



# MY OASIS MARKETING MATERIALS

OASIS provides a broad mix of informational materials in the form of brochures, posters, and videos for use in your practice and to distribute to your consumers/patients as new members join your community.

	Product Name	Ref #
	My OASIS How to Register on a Smart Phone/Tablet Booklet	LIT-MOL-RegM
	My OASIS How to Register on a Desktop Computer Booklet	LIT-MOL-RegD
	My OASIS Member Registration Postcard	LIT-MOL-PC
	My OASIS Member Registration Wall Chart	LIT-MOL-WC
	My OASIS Doctor Recommendation Pad	LIT-OTC-MORX





# LUBRICANT EYE DROPS



## RELIEF FOR DRY EYES

Oasis TEARS Lubricant Eye Drops' proprietary formulation soothes, lubricates, and relieves burning, irritation, dryness, and other symptoms related to dry eye. Oasis TEARS® line of dry eye relief products are available through eye care professionals and not available in major retail stores. Oasis TEARS is an over-the-counter lubricant eye drop.

## DETAILS

- Dry eye relief
- Moistening, lubrication and comfort
- For use as needed
- 3 treatment options
- Single use containers or multi-use bottle

Ask how Oasis can support your practice with dry eye solutions implementation and staff training. **NOT AVAILABLE FOR EXPORT.**

Product Name	Description	Qty Per Carton	Ref #
Oasis TEARS PF Plus	Preservative-free, 0.34 fl. oz. (10ml) bottle	20 bottles	OT6500
Oasis TEARS PF	Preservative-free, 0.34 fl. oz. (10ml) bottle	20 bottles	OT6400
Oasis TEARS	Preservative-free, 30 single-use containers	45 boxes	OT6100
Oasis TEARS PLUS	Preservative-free, 30 single-use containers	45 boxes	OT6200
Oasis TEARS	With a mild preservative, 0.34 fl. oz. (10ml) bottle	20 bottles	OT6310





# HOT & COLD EYE MASK

## RELIEF FOR DRY EYES

Oasis REST & RELIEF provides heat therapy for the relief of dry eye symptoms. Use as directed.



## PRODUCT FEATURES

- Temperature retaining beads that conform to the body in a soft knit cloth mask for maximum comfort
- Adjustable strap for comfortable compression
- Washable cover

Product Name	Description	Qty Per Carton	Ref #
Oasis REST & RELIEF EYE MASK	Eye mask for hot and cold use, 1 mask and 1 mask cover per box	24 Masks	OM8000
Oasis REST & RELIEF EYE MASK Cover	Eye mask cover, 1 mask cover	1 Mask Cover	OM8001

"My ophthalmologist suggested that I try Oasis Tears Eye Drops a few weeks ago. I have been having a lot of trouble with dry eyes & consequent diminished vision since I had Lasik eye surgery about 3 years ago. I have tried several different brands of eye drops without getting relief that I need. The Oasis drops were different right away. It was amazing to only need to use eye drops 4 times a day instead of 6-10 times a day. They have relieved the dryness in my eyes much better than other drops and my vision has improved in the past couple of weeks. I also am very pleased that these drops are small so I don't have tears running down my cheeks and a recovery time before I can see. These drops are so easy, work better and are much more convenient than others. I have also recommended these drops to family members who have trouble with dry eyes. So glad to be able to use this better product. Wish they were available in stores also. Thank you."

**Glenna C., California**

"I love the tears! I really think your tear is one of the best, if not the best on the market! I am following your suggestion of placing a drop on the patient's eye while in the exam chair. They can tell a real difference. They really do last long on the eye."

**Michael H., California**

"...I especially love the eye mask! I use it every day. It also helps before or after my Restasis drops, because the drops burn. I highly recommend the mask."

**Norma H., California**





Oasis LID & LASH®  
+ tea tree oil  
CLEANSE + HYDRATE

## EYELID & LASH CLEANSING PADS

### DAILY LID & LASH HYGIENE

Oasis TEARS Oasis LID & LASH + tea tree oil cleanses and hydrates eyelids and lashes with our new gentle tea tree formulation to maintain eyelid and lash appearance. Our soft pads gently remove eye debris and dust particles to cleanse delicate skin around your eyes. Leave on to air-dry smooth for a refreshing experience of clean with tea tree.



NEW!

### DETAILS

- Gentle for daily use
- Paraben Free
- 60 cleansing pads

Product Name	Description	Qty Per Carton	Ref #
Oasis TEARS Oasis LID & LASH + tea tree oil	Eyelid and lash cleansing pad, 60 pre-soaked pads	24 jars	OT7300



HYPOCHLOROUS®  
Eyelid Cleansing Spray

## EYELID CLEANSING SPRAY

### DAILY LID & LASH HYGIENE

Oasis TEARS HYPOCHLOROUS is a gentle, non-toxic and natural antimicrobial formula for use daily on eyelids and lashes. Oasis TEARS HYPOCHLOROUS contains Hypochlorous Acid a known antimicrobial, which serves to inhibit the growth of microorganisms. Hypochlorous Acid cleanses irritated eyelids of bioburden. This natural antimicrobial can be sprayed directly on closed eyelids to soothe them from irritation. The no-rinse spray can air-dry on to manage bioburden.



NEW!

### DETAILS

- Gentle for daily use
- .012% Hypochlorous Acid
- Natural Antimicrobial

Product Name	Description	Qty Per Carton	Ref #
Oasis TEARS HYPOCHLOROUS Spray	Eyelid and lash cleansing spray, 1.86 fl. oz. (55ml) bottle	20 bottles	OT7400







## DIETARY SUPPLEMENT

### Oasis TEARS OMEGA<sup>®</sup> 3

Oasis TEARS OMEGA 3 uses proprietary methods to remove toxins and metals for highest purity, plus an additional process to enhance the concentration of EPA and DHA.

#### DETAILS

- Tryglyceride form
- Total DHA & EPA 1000mg per softgel
  - DHA 429mg per softgel\*\*
  - EPA 571mg per softgel\*\*
- 90 softgels = 3 softgels per serving
- Natural lemon flavor added
- No fish taste

\*\*Total including other omegas per softgel = 1143mg

Product Name	Description	Qty Per Carton	Ref #
Oasis TEARS OMEGA3	Dietary Supplement	12 Bottles	<b>ON3000</b>



LID HYGIENE & DIETARY SUPPLEMENTS



## DIETARY SUPPLEMENT

### Oasis TEARS VISION<sup>®</sup>

Oasis TEARS VISION's patent pending formulation supports dietary supplementation for general eye health.

#### DETAILS

- Patent pending formulation
  - Maqui Berry (30mg)
  - Astaxanthin (6mg)
  - Zeaxanthin (5mg)
- Lutein (25mg)
- DHA (200mg)
- Small softgel capsule
- 30 softgels = 1 softgel per serving

Product Name	Description	Qty Per Carton	Ref #
Oasis TEARS VISION	Dietary Supplement	12 Bottles	<b>ON3010</b>



Oasis TEARS OMEGA 3 and Oasis TEARS VISION are free from \*gluten, \*dairy, and made with NON-GMO ingredients. Oasis TEARS OMEGA 3 and Oasis TEARS VISION are produced in a NSF, FDA Registered, and Inspected Facility. \*Oasis TEARS OMEGA 3 and Oasis TEARS VISION are produced in a facility that processes other ingredients that may contain these allergens.



# DRY EYE RELIEF KITS

## OASIS DRY EYE RELIEF KIT 1: K-DEKITS-003

Product Name	Description	Dimensions	Qty Per Bag
Oasis TEARS PF Plus	Preservative-free eye drops		1 bottle
Oasis TEARS Oasis LID & LASH + tea tree oil	Eyelid and lash cleansing pad, 60 pre-soaked pads		1 jar
Oasis REST & RELIEF EYE MASK	Eye mask for hot and cold use		1 mask and 1 mask cover
OASIS Dry Eye Kitting Bag	Fabric Bag with Zipper	9 x 3 x 6.5	1 bag



## OASIS DRY EYE RELIEF KIT 2: K-DEKITS-004

Product Name	Description	Dimensions	Qty Per Bag
Oasis TEARS PF Plus	Preservative-free eye drops		1 bottle
Oasis TEARS Oasis LID & LASH + tea tree oil	Eyelid and lash cleansing pad, 60 pre-soaked pads		1 jar
Oasis REST & RELIEF EYE MASK	Eye mask for hot and cold use		1 mask and 1 mask cover
Oasis TEARS OMEGA 3	Dietary Supplement, 90 Softgels		1 bottle
OASIS Dry Eye Kitting Bag	Fabric Bag with Zipper	9 x 3 x 6.5	1 bag



# DRY EYE RELIEF KITS

## OASIS DRY EYE RELIEF KIT 3: K-DEKITS-005

Product Name	Description	Dimensions	Qty Per Bag
Oasis TEARS PF	Preservative-free eye drops		1 bottle
Oasis TEARS HYPOCHLOROUS Spray	Eyelid and lash cleansing spray		1 bottle
Oasis REST & RELIEF EYE MASK	Eye mask for hot and cold use		1 mask and 1 mask cover
OASIS Dry Eye Kitting Bag	Fabric Bag with Zipper	9 x 3 x 6.5	1 bag



DRY EYE RELIEF KITS

## CONTACT US TO CREATE A CUSTOM VISION CARE KIT

Pick from our OTC line of dry eye solutions to create your practice's custom pre-op and post-op kits using a choice of canvas bag or box with insert. Dimensions are approximate.

Product Name	Description	Dimensions	Ref #
OASIS Dry Eye Kitting Bag	Fabric Bag with Zipper	9 x 3 x 6.5	S7100
Oasis TEARS Box	Box without Insert	11 x 8 x 3	S8100



# FORM FIT<sup>®</sup>

## HYDROGEL CANALICULAR PLUGS

### LONG TERM OCCLUSION





FORM FIT is truly a one size plug solution. Once inserted into the canaliculus, FORM FIT hydrates in approximately 10 minutes. As it hydrates, FORM FIT increases in size until it completely fills the vertical canalicular cavity.

### HOW IT WORKS

FORM FIT is made of a safe and biocompatible hydrogel material that expands into a soft, pliable, gelatinous material when it contacts tear film. FORM FIT securely conforms to the vertical canaliculus so there is no need for a cap on the end of the plug. This greatly increases patient comfort and decreases foreign body sensation.

### DETAILS

- Hydrates in approximately 10 minutes
- Easily removed via irrigation of the tear duct
- Biocompatible, no foreign body sensation
- Sterile
- Conveniently pre-loaded in inserter
- 0.3mm, one size fits all



Product Name	Qty Per Box	Diameter	Ref #
FORM FIT Hydrogel Canalicular Plug	10pk	0.3mm	6303
FORM FIT Hydrogel Canalicular Plug	2pk	0.3mm	6303-T






# SILICONE PUNCTUM PLUGS

## LONG TERM OCCLUSION

The distinctive low-profile cap combined with a softer, flexible silicone maximizes patient comfort. Now available with two inserter options. SOFT PLUG Silicone Punctum Plugs are removable with forceps per instruction for use.



## DETAILS

- **NEW** Pen-Style Inserter with dilator tip - 2pks, 6pks
- Original Inserter with dilator tip - 2pks, 6pks, 20pks
- Sterile, preloaded on inserter in 2-pack or 6-pack dispenser box
- Low-profile cap intended to reduce dislocation
- Economy pack includes 20 plugs and two non-sterile inserters. Holders equipped with an alignment slot for quick loading



Product Name	Qty Per Box	Diameter	Original Inserter Ref #	Pen-Style Inserter Ref #
SOFT PLUG Silicone, Preloaded, Micro	2pk	0.4mm	6610	-
SOFT PLUG Silicone, Preloaded, Mini	2pk	0.5mm	6611	-
SOFT PLUG Silicone, Preloaded, Petite	2pk	0.6mm	6612	-
SOFT PLUG Silicone, Preloaded, Small	2pk	0.7mm	6613	-
SOFT PLUG Silicone, Preloaded, Medium	2pk	0.8mm	6614	-

Product Name	Qty Per Box	Diameter	Original Inserter Ref #	Pen-Style Inserter Ref #
SOFT PLUG Silicone, Preloaded, Micro	6pk	0.4mm	6610-D	NEW! 7710-D
SOFT PLUG Silicone, Preloaded, Mini	6pk	0.5mm	6611-D	NEW! 7711-D
SOFT PLUG Silicone, Preloaded, Petite	6pk	0.6mm	6612-D	NEW! 7712-D
SOFT PLUG Silicone, Preloaded, Small	6pk	0.7mm	6613-D	NEW! 7713-D
SOFT PLUG Silicone, Preloaded, Medium	6pk	0.8mm	6614-D	NEW! 7714-D

Product Name	Qty Per Box	Diameter	Original Inserter Ref #
SOFT PLUG Silicone, Economy Pack, Micro	20pk	0.4mm	6610-EP
SOFT PLUG Silicone, Economy Pack, Mini	20pk	0.5mm	6611-EP
SOFT PLUG Silicone, Economy Pack, Petite	20pk	0.6mm	6612-EP
SOFT PLUG Silicone, Economy Pack, Small	20pk	0.7mm	6613-EP
SOFT PLUG Silicone, Economy Pack, Medium	20pk	0.8mm	6614-EP






## SILICONE PUNCTUM PLUG, FLOW CONTROL

### LONG TERM OCCLUSION WITH PARTIAL DRAINAGE

The distinctive low-profile cap combined with a softer, flexible silicone maximizes patient comfort. The inserter minimizes potential for accidental release of the plug. Plug is designed with a center channel allowing for partial tear flow. SOFT PLUG Flow Control Plugs are removable with forceps (per instruction for use).



### DETAILS

- Available in 1pk and 2pks
- Silicone
- Sterile, preloaded on inserter
- For partial occlusion with limited tear drainage


Product Name	Qty Per Box	Diameter	Ref #
SOFT PLUG Flow Control, Preloaded, Petite	2pk	0.6mm	6612-FC
SOFT PLUG Flow Control, Preloaded, Small	2pk	0.7mm	6613-FC
SOFT PLUG Flow Control, Preloaded, Medium	2pk	0.8mm	6614-FC



## COLLAGEN INTRACANALICULAR PLUGS

### TEMPORARY OCCLUSION

The Collagen Plugs have a smooth proprietary finish that simplifies insertion of the plug. Collagen Plugs are ideal for short-term, diagnostic, and post-surgical occlusion.



### DETAILS

- Degradable material lasts 2-5 days
- Proprietary finish simplifies insertion
- Evaluate patient acceptance
- 2mm length
- Sterile
- Insert with forcep

Product Name	Qty Per Box	Diameter	Ref #
SOFT PLUG Collagen Intracanalicular Plug	60pk	0.2mm	6502
SOFT PLUG Collagen Intracanalicular Plug	60pk	0.3mm	6503
SOFT PLUG Collagen Intracanalicular Plug	60pk	0.4mm	6504






# EXTENDED DURATION PLUGS


## MEDIUM TERM OCCLUSION FOR UP TO 3 OR 6 MONTHS

SOFT PLUG Extended Duration Punctum Plugs are degradable canalicular occluders that are an effective option for transient dry eye symptoms before or after surgery, seasonally, etc. and for patients seeking to experience non-permanent options of punctal occlusion.

**SOFT PLUG® EXTENDED DURATION  
180 CANALICULAR PLUGS**  
Up to 180 days.



**SOFT PLUG® EXTENDED  
DURATION PLUGS**  
Up to 90 days.



### DETAILS

- Pre and post-surgical use
- Sterile
- Insert with forceps
- 2mm length

Product Name	Qty Per Box	Diameter	Ref #
SOFT PLUG Extended Duration 180 Canalicular Plug	20pk	0.3mm	<b>6203-A</b>
SOFT PLUG Extended Duration 180 Canalicular Plug	20pk	0.4mm	<b>6204-A</b>
SOFT PLUG Extended Duration 180 Canalicular Plug	20pk	0.5mm	<b>6205-A</b>
SOFT PLUG Extended Duration 180 Canalicular Plug	40pk	0.3mm	<b>6203</b>
SOFT PLUG Extended Duration 180 Canalicular Plug	40pk	0.4mm	<b>6204</b>
SOFT PLUG Extended Duration 180 Canalicular Plug	40pk	0.5mm	<b>6205</b>
SOFT PLUG Extended Duration 180 Canalicular Plug, Variety Pack	20pk	3 pair 0.3mm, 4 pair 0.4mm, 3 pair 0.5mm	<b>6200V</b>
SOFT PLUG Extended Duration Plug	20pk	0.2mm	<b>6402-A</b>
SOFT PLUG Extended Duration Plug	20pk	0.3mm	<b>6403-A</b>
SOFT PLUG Extended Duration Plug	20pk	0.4mm	<b>6404-A</b>
SOFT PLUG Extended Duration Plug	20pk	0.5mm	<b>6405-A</b>
SOFT PLUG Extended Duration Plug	40pk	0.2mm	<b>6402</b>
SOFT PLUG Extended Duration Plug	40pk	0.3mm	<b>6403</b>
SOFT PLUG Extended Duration Plug	40pk	0.4mm	<b>6404</b>
SOFT PLUG Extended Duration Plug	40pk	0.5mm	<b>6405</b>
SOFT PLUG Extended Duration Plug, Variety Pack	20pk	1 pair 0.2mm, 3 pair 0.3mm, 4 pair 0.4mm, 2 pair 0.5mm	<b>6400V</b>




# ACCESSORIES & DRY EYE DIAGNOSTICS



S6100



6648



4003-IK

ACCESSORIES & DRY EYE DIAGNOSTICS

Product Name	Description	Qty Per Box	Ref #
Tear Flo Test Strips	Sterile strip	100	6512

Product Name	Description	Qty Per Box	Ref #
SOFT PLUG Punctal Plug Gauge	0.5mm and 0.6mm tips	1	6625
SOFT PLUG Punctal Plug Gauge	0.7mm and 0.8mm tips	1	6627
SOFT PLUG Punctal Plug Gauge Set	Includes both instruments and sterilization tray	1 set	6630

Product Name	Description	Qty Per Box	Ref #
Straight Jewelers Forceps	Straight #3	1	6643
Curved Jewelers Forceps	Curved #7	1	6647
Grooved Forceps	For the insertion of Extended Duration & Collagen Plugs	1	6648
Castroviejo Lacrimal Dilator	Double End	1	6640

Product Name	Description	Qty Per Box	Ref #
Lacrimal Irrigation Kit	One 3mL syringe and one 25g x 1 1/8" lacrimal cannula, curved, open end	10	4003-IK
Lacrimal Intubation Set	23g	3	6001
Lacrimal Intubation Set with Retrieval Device	23g	3	6004
Lacrimal Intubation Kit	27g	3	6004-27
Oasis TEARS Bag	Clear, plastic bag with Oasis TEARS logo	100 bags	S6100















# MARKETING MATERIALS

OASIS® offers educational materials in the form of brochures, displays, posters, and videos for use in your practice and to distribute to your patients.

To view and download digital versions of all resources, go to Practice & Patient Materials at [www.oasismedical.com](http://www.oasismedical.com)

Product Name	Ref #	Product Name	Ref #
	Dry Eye Catalog <b>LIT-USDE-CAT</b>		Oasis TEARS Doctor Recommendation Pad, with OASIS Contact Information <b>LIT-OTRX-W</b>
	Oasis TEARS Brochure <b>6544C</b>		Oasis TEARS Doctor Recommendation Pad, without OASIS Contact Information <b>LIT-OTRX-NB</b>
	Punctal Plug Brochure <b>LIT-PP-PB</b>		Dry Eye Symptoms Wall Chart <b>LIT-DES-WC</b>
	Dry Eye Products Patient Brochure <b>LIT-DEP-PB (English)</b> <b>LIT-DEP-PBS (Spanish)</b> <b>LIT-DEP-PBC (Trad. Chinese)</b>		Punctal Plug Wall Chart <b>LIT-WCDE</b>
	Oasis TEARS Point-of-Purchase Display <b>LIT-OT-POP</b>		Surgical Catalog <b>LIT-SRG-CAT</b>

MARKETING MATERIALS



# ORDERING INFORMATION

Orders can be placed 24 hours via our website, fax, or the phone message center.

## CUSTOMER SERVICE HOURS

7:00 AM - 4:30 PM (PST),

Monday - Friday, Excluding holidays.

**USA:** 844.820.8940

800.631.7210 FAX

**INTERNATIONAL:** 909.305.5400

909.305.9987 FAX

**PAYMENT TERMS:** NET 30 DAYS

**SHIPPING:** FOB GLENDORA, CALIFORNIA

Insurance and shipping charges added to invoice.

international@oasismedical.com

**WEB (US ONLY):** www.oasismedical.com

## RETURN MERCHANDISE AUTHORIZATION (RMA)

Unopened items returned for credit must be accompanied by an RGA (return goods authorization number) and the reason for the return. An RMA can be requested through customer service.

Allowance Schedule for acceptable returns	
Sellable merchandise returned in factory condition within 60 days	Full Credit
Sellable merchandise returned in factory condition within 60-90 days	Cost less 40%
Sellable merchandise returned in factory condition 90 days to 1 year	Cost less 60%
Merchandise returned one year after invoice date	No Credit
Merchandise unsuitable for resale	No Credit
Merchandise defective at time of receipt by customer	Full Credit
Merchandise shipped in error	Full Credit

### Merchandise that will not be accepted:

- “Special Order” or “Custom Products”
- Merchandise not in standard OASIS® packaging
- Opened products
- Discontinued products

*NOTE: Customer must prepay all freight charges on stock returns except for defective material or merchandise shipped in error.*

## DOMESTIC SHIPPING POLICY

Shipments from OASIS® are FedEx only. Choices are Ground, 3-Day Select, 2-Day, Priority Overnight (by 10:30am) or Standard Overnight (by 4:30pm). Shipping costs will be added when order is placed. Customer’s FedEx or UPS account number may be used if provided to Customer Service at time order is placed. We cannot ship to PO boxes. We cannot be responsible for delays in transit. In-stock product orders received before 2pm PST Monday-Friday, excluding holidays, will usually ship on the same business day. Prepayment orders will ship once payment is received. Custom orders will require additional shipment times and will ship when the product is available.

## WEBSITE ORDERS

Visit [www.oasismedical.com](http://www.oasismedical.com) to see terms of use, privacy policy, and returns authorization by web.















# STERILIZATION METHOD

Products are sterilized by ethylene oxide, gamma radiation, electron beam, or steam according to EN 556 (European) guidelines and ANSI/AAMI/ISO 11135, 11137 or 17665 standards for medical device sterilization as applicable, to a sterility assurance level (SAL) of 10<sup>-6</sup>.








## PRODUCT LABELING

OASIS® uses internationally recognized symbols on its product labeling.

	European Union CE Mark and notified body ID number
	Consult instructions for use.
<b>Rx Label</b>	Physician use only
	Do not reuse, single use
<b>LOT</b> LSmmyyA	Lot or batch code
<b>MD</b>	Medical Device
 yyyy-mm-dd	Use by. Placed next to Expiration Date – year/month
<b>REF</b> [Part No.]	Reference Catalog Number (Numbers/Letter represent part number i.e. 6612D)
<b>EC</b>   <b>REP</b>	European Authorized Representative
	Manufacturer
<b>UDI</b>	UDI
	Patient Information Website
	Contains biological material of animal origin
<b>CH</b>   <b>REP</b>	Swiss Authorized Representative

	Do not resterilize
	Do no use if package is damaged
	No Latex
	Keep dry
	Temperature limit

## METHODS OF STERILIZATION

	Gamma or E-beam Radiation, Single Sterile Barrier System
	Gamma or E-beam Radiation, Double Sterile Barrier System
	Gamma or E-beam Radiation, Single Sterile Barrier System with Protective Packaging Inside
	Gamma or E-beam Radiation, Single Sterile Barrier System with Protective Packaging Outside
	Ethylene Oxide, Single Sterile Barrier System
	Ethylene Oxide, Double Sterile Barrier System
	Steam or Dry Heat, Single Sterile Barrier System



"I have Sjögren's and dry eyes are one of my symptoms. I literally had dry patches on my eyes. I needed to have a product I can depend on to keep them hydrated and stop that gritty itchy feeling. I know it may sound cliché, but it's true! My eyes are the desert and you are my Oasis! I've tried so many products, which means I've spent a lot of money. Not anymore. I've found a product that not only feels good when I put it in, but continues to feel good. I was so frustrated, but not anymore. Thank you OASIS!"

**Charlene D., New Jersey**

---

"I love your product. It is the only thing that seems to help my eyes since my surgery 3 months ago. After my epi-lasik I had numerous erosions and was in total discomfort no matter what the eye doctor gave me for eye relief. I asked if he didn't have any other new products that I could try that might give me more relief for the dryness in my eyes. He said he had a new product... Oasis TEARS PLUS. From the first day I used it my eyes have felt so much better. I am at the point one month later that I don't even feel the dryness everyday when I wake up. Thanks for your product."

**Patti M., Minnesota**

Ask your OASIS® representative about our Dry Eye Solutions.



Glendora, CA, USA  
(844) 820-8940 (USA Toll Free) | (909) 305-5400  
oasis@oasismedical.com | www.oasismedical.com

LIT-USDE-CAT Rev14 3/2022